

Troubleshooting Guide for Transfection




Table 11.6. Troubleshooting Guide for Transfection.

Problem	Possible cause and recommended solution
1. Low transfection efficiency	<p>1.1. Inefficient complex formation. Always vortex the mixture immediately after the addition of the reagent to DNA.</p> <p>1.2. Suboptimal reagent/DNA ratio. Optimize the quantity of transfection reagent added to the fixed amount of DNA.</p> <p>1.3. Suboptimal quantity of DNA. Optimize the amount of DNA used for transfection. Keep the transfection reagent/DNA ratio constant</p> <p>1.4. Poor polyplex/cell surface contact. Centrifuge the culture plates, if it will not harm the cells. Caution – centrifugation may harm some primary cells.</p> <p>1.5. Poor DNA quality. Use high quality DNA with an A_{260}/A_{280} ratio greater than 1.8.</p> <p>1.6. Suboptimal cell confluency. Optimize cell plating conditions. Ensure that adhered cells are 50-70% confluent at the time of transfection. Ensure that suspension cells are in logarithmic growth phase at the time of transfection.</p> <p>1.7. Mycoplasma contamination. Mycoplasma infection in cell culture often results in poor and/or non-reproducible transfection. Regularly check your cells for mycoplasma infection.</p>
2. High cellular toxicity	<p>2.1. Toxic transgene. Verify if the expressed transgene is toxic</p> <p>2.2. Suboptimal incubation conditions. Reduce incubation time of the polyplexes with the cells. Replace the transfection mixture 3-6 hours later with fresh growth medium.</p> <p>2.3. Suboptimal quantity of DNA. Reduce the quantity of DNA used for transfection</p> <p>2.4. Cell density is too low. Increase the plating density of cells used for transfection.</p>