




## Troubleshooting Guide



**Table 1.37.** Troubleshooting guide for DNA digestion.

Problem	Possible cause and recommended solution
<p><b>1. Incomplete digestion or no digestion</b></p>	<p><b>Assess enzyme activity</b> The restriction enzyme may lose activity due to improper storage or handling. Perform a digestion reaction with 1 µg of a standard control DNA, e.g. Lambda DNA (<i>dam</i><sup>-</sup>, <i>dcm</i><sup>-</sup>) (#SD0021).</p> <p><b>1.1. Inactive enzyme.</b> If the enzyme does not cut the control DNA:</p> <ul style="list-style-type: none"> <li>• Check the expiration date.</li> <li>• Verify that the enzyme has been stored at -20°C.</li> <li>• Check the temperature of your freezer. Do not allow the temperature go below -20°C as the enzyme may freeze and multiple freeze thaw cycles (more than 3 cycles) may result in reduced enzyme activity.</li> </ul> <p><b>1.2. Suboptimal reaction conditions.</b></p> <p><b>1.2.1. Suboptimal digestion protocol.</b></p> <ul style="list-style-type: none"> <li>• Follow digestion protocol specified for the restriction enzyme and type of substrate DNA. All FastDigest® enzymes are experimentally tested on Lambda DNA (or other control substrate DNA), plasmid and genomic DNA as well as PCR products. Please refer to Table 1.3 on p.71 for specific recommendations.</li> <li>• For FastDigest® enzymes use FastDigest® or FastDigest® Green Buffer. All FastDigest® enzymes are 100% active in these buffers.</li> <li>• For conventional restriction enzymes use the recommended reaction buffer supplied with the enzyme. For double digestions follow the recommendations of the DoubleDigest™ engine at <a href="http://www.fermentas.com/doubledigest">www.fermentas.com/doubledigest</a>.</li> <li>• Use additives where required.</li> <li>• Perform the reaction at the optimal temperature specified for the restriction enzyme. For double digestions with enzymes requiring different incubation temperatures perform sequential DNA cleavage: complete the first digestion reaction at the lower temperature, add the second enzyme and increase the digestion temperature for the second enzyme cleavage.</li> <li>• Ensure the volume of the reaction mixture was not reduced due to evaporation during incubation; the increase in salt concentration may reduce enzyme activity. For thermophilic enzymes use a heat block with a hot bonnet, e.g. a PCR cycler.</li> </ul> <p><b>1.2.2. Improper dilution of conventional enzyme.</b></p> <ul style="list-style-type: none"> <li>• Dilute conventional restriction enzymes with Dilution Buffer for Restriction Enzymes (#B19). Restriction enzymes diluted with this buffer are stable for at least 3-4 weeks at -20°C.</li> <li>• Never dilute enzymes in water or 10X reaction buffer.</li> <li>• Never dilute enzymes in 1X reaction buffer in the absence of DNA.</li> <li>• Use the recommended amount of FastDigest® enzymes, which are experimentally tested on four different DNA substrates (<i>see</i> Table 1.3 on p.71).</li> </ul> <p><b>1.2.3. Improper reaction assembly.</b></p> <ul style="list-style-type: none"> <li>• The restriction enzyme should always be the last component added to the reaction mixture.</li> <li>• The restriction enzyme may be inactivated if added directly to a 10X reaction buffer.</li> </ul> <p><b>1.2.4. Excess glycerol in the reaction mixture.</b></p> <ul style="list-style-type: none"> <li>• The glycerol concentration in the reaction mixture should not exceed 5%. Thus, the volume of the restriction enzyme added to the mixture should not exceed 1/10 of the total reaction volume.</li> <li>• Enzymes sensitive to high glycerol concentration include: Alw21I, Bpil, Bsp68I, BspT1, Eco32I, Eco91I, EcoRI, Hin6I, HinfI, Mph1103I, Mva1269I and NcoI.</li> </ul> <p><b>1.2.5. Suboptimal DNA concentration.</b> The optimal range of DNA concentration in the reaction mixture is 0.02-0.1 µg/µl.</p> <p><b>1.3. Unsuitable DNA template or contaminated DNA solution.</b> If the enzyme is active in the control digest, assay the substrate DNA solution for inhibitory contaminants in a mixing experiment with control template, e.g. Lambda DNA (<i>dam</i><sup>-</sup>, <i>dcm</i><sup>-</sup>) (#SD0021). Perform a control digest with two templates, control and sample, in one reaction mixture. Do not exceed the optimal DNA concentration in the reaction mixture (0.02-0.1 µg/µl).</p> <ul style="list-style-type: none"> <li>• The sample template is contaminated if neither the control DNA nor sample DNA template is digested (<i>see</i> 1.3.1).</li> <li>• The sample template is not contaminated if the control DNA template is digested but the sample template is not. Poor digestion of the experimental template is caused by errors in the DNA sequence (<i>see</i> 1.3.2), methylation effects (<i>see</i> 1.3.3) or structure of the DNA substrate (<i>see</i> 1.3.4).</li> </ul> <p><b>Note</b> Always ensure that the control DNA contains a recognition site for the enzyme present in the reaction. For example, there is no NotI recognition site in lambda DNA.</p>

(continued on next page)




**Table 1.37.** Troubleshooting guide for DNA digestion.

Problem	Possible cause and recommended solution
1. Incomplete digestion or no digestion	<p><b>1.3.1. Contaminants in the DNA solution.</b></p> <ul style="list-style-type: none"> <li>• Template DNA may contain residual SDS, EDTA, proteins, salts or nucleases. Repurify the template using a GeneJET™ PCR Purification Kit (#K0701) or by phenol/chloroform extraction and ethanol precipitation (see p.358). DNA A<sub>260/280</sub> ratio should be 1.8-2.0. To remove EDTA and salts, wash the pellet with 70% cold ethanol.</li> <li>• For reliable and reproducible plasmid miniprep purity, use the GeneJET™ Plasmid Miniprep Kit (#K0502).</li> <li>• For digestion of unpurified PCR products, dilute DNA at least 3-fold in the recommended 1X restriction enzyme buffer see protocols on p.70 or p.160.</li> <li>• If the template DNA has been purified using silica or resin suspensions, remove all remaining particles by centrifugation for 10 min at 10,000 rpm and ensure that no resin is carried over while transferring the DNA solution into a new tube.</li> </ul> <p><b>1.3.2. The substrate DNA does not contain a recognition sequence for the restriction enzyme.</b></p> <ul style="list-style-type: none"> <li>• Re-check the DNA sequence and cloning strategy.</li> <li>• Determine if the restriction enzyme selected requires more than one site per target DNA for 100% activity (see also 1.3.4.3).</li> <li>• Check literature for known site preferences for the restriction enzyme (see also 1.3.4.4).</li> <li>• If the recognition sequence was introduced by PCR primers, verify that the primer sequence contains the recognition site.</li> </ul> <p><b>1.3.3. Methylation effects.</b></p> <p>Restriction enzyme may be inhibited by methylation of the recognition site.</p> <ul style="list-style-type: none"> <li>• Identify which type of DNA methylation can occur on the recognition site and determine if the methylation impairs or blocks DNA digestion with the enzyme. See Digestion of Methylated DNA on p.175 and use Tables 1.13 and 1.20 on pp.176-181.</li> </ul> <p>If methylation impairs or blocks DNA cleavage:</p> <ul style="list-style-type: none"> <li>• Propagate your plasmid in an <i>E.coli dam<sup>-</sup>, dcm<sup>-</sup></i> strain (the <i>E.coli</i> GM2163 <i>dam<sup>-</sup>, dcm<sup>-</sup></i> strain; #M0099, is available upon request with the purchase of any Fermentas product),</li> <li>• Use the REsearch™ engine at <a href="http://www.fermentas.com/research">www.fermentas.com/research</a> or check the Fermentas catalog for the availability of a restriction enzyme isoschizomer not sensitive to DNA methylation.</li> <li>• Using of restriction enzyme which requires a methylated recognition sequence (DpnI or SgeI) for digestion of unmethylated DNA will result in no DNA cleavage. Propagate your plasmid in <i>E.coli dam<sup>+</sup> or dcm<sup>+</sup></i> strains (please refer to p.486 for genotype information of some common <i>E.coli</i> strains) to get DNA with methylated DpnI or SgeI recognition sequences. Alternatively, in the case of DpnI, the neoschizomers Bsp143I or MboI can be used to digest non-methylated DpnI recognition sites.</li> </ul> <p><b>Note</b></p> <p>When PCR is carried out with standard dNTPs and non-methylated primers the resulting DNA product is not methylated.</p> <p><b>1.3.4. Structure of substrate DNA.</b></p> <p><b>1.3.4.1. Supercoiled plasmid DNA.</b></p> <ul style="list-style-type: none"> <li>• Use FastDigest® enzymes which are qualified for supercoiled DNA and provided with specific recommendations for each enzyme (see Table 1.3 on p.71).</li> <li>• For some conventional restriction enzymes, additional units are required to completely digest supercoiled plasmids (e.g. 5-10 u (1 µl) of restriction enzyme per 1 µg of DNA). Check notes in the catalog description of the enzyme or refer to the Certificate of Analysis.</li> </ul> <p><b>1.3.4.2. Proximity of the recognition sequence to the DNA ends.</b></p> <p>Some restriction enzymes cleave DNA poorly if the recognition site is too close to the end of the DNA molecule.</p> <ul style="list-style-type: none"> <li>• For FastDigest® enzymes refer to Table 1.3 on p.71 or the product description to determine the effectiveness of restriction enzyme cleavage at the ends of DNA.</li> <li>• For conventional restriction enzymes refer to Tables 1.10 (p.171) and 1.11 (p.172).</li> </ul> <p><b>1.3.4.3. Restriction enzyme requires at least two sites per DNA molecule for optimal activity.</b></p> <p>Some restriction enzymes such as AarI, BvuI, Cfr10I, Cfr42I, Eco57I, EcoRII, LweI, SfiI require at least two target sites per DNA molecule for efficient cleavage (for more details see Site Preferences by Restriction Enzymes on p.173). If there is only one recognition site per DNA molecule, add a DNA oligonucleotide containing the recognition site.</p> <p><b>1.3.4.4. Site Preferences by Restriction Enzymes.</b></p> <p>The DNA sequence surrounding the recognition site may influence the efficiency of digestion. Some DNA sites are cleaved slowly or not cleaved at all (for more details see p.173) due to the surrounding sequence. Use additional units (5-10 u) of conventional restriction enzyme per 1 µg of DNA or determine if an isoschizomer has superior cleavage efficiency (see Table 1.25 on p.207 or REsearch™ engine at <a href="http://www.fermentas.com/research">www.fermentas.com/research</a>).</p>

(continued on next page)

**Table 1.37.** Troubleshooting guide for DNA digestion.

Problem	Possible cause and recommended solution
1. Incomplete digestion or no digestion	<p><b>1.4. Water contains impurities.</b></p> <ul style="list-style-type: none"> <li>Compare your results using commercially available nuclease free, molecular biology grade water, e.g. Water, nuclease-free (#R0581). Check the quality of the water used in your lab.</li> <li>Check the pH and conductivity of water. The pH of high quality water should be 5.5-6.0 with a resistance of <math>\geq 18 \text{ M}\Omega</math>.</li> <li>Centrifuge (10 min, 10,000 rpm) 1 ml of water and check if there is a visible pellet.</li> <li>Determine if water contains nucleases or bacterial contamination (see 3.2 for control reactions).</li> </ul>
2. Unexpected cleavage pattern	<p><b>2.1. Star activity (relaxed specificity) of restriction enzyme</b> (see p.174 for more details).</p> <ul style="list-style-type: none"> <li>Use FastDigest® restriction enzymes. For these enzymes the incubation time without star activity is determined.</li> <li>Reduce the amount of conventional restriction enzyme (use no more than 10 u).</li> <li>Use the recommended reaction buffer.</li> <li>Ensure that the glycerol concentration in the reaction mixture does not exceed 5%.</li> <li>Reduce the incubation time. For FastDigest® enzymes – refer to Table 1.3 on p.71 for maximum incubation times.</li> <li>Ensure the volume of the reaction mixture was not reduced due to evaporation during incubation. The resulting increase in glycerol concentration may cause star activity.</li> </ul> <p><b>2.2. Contamination with another restriction enzyme.</b> The restriction enzyme or buffer may be contaminated with another restriction enzyme due to improper handling. Use a new tube of enzyme and/or buffer.</p> <p><b>2.3. Contamination with another substrate DNA.</b> The sample DNA may contain a mixture of two or more different DNAs. Prepare new sample of DNA.</p> <ul style="list-style-type: none"> <li>For plasmid DNA preparation pick one isolated colony of recombinant <i>E.coli</i> grow cells and purify plasmid with GeneJET™ Plasmid Miniprep Kit (#K0503).</li> <li>For PCR products: check the product purity on an agarose gel. If necessary, purify the PCR product prior to digestion with GeneJET™ Gel Extraction Kit (#K0692) or Silica Bead DNA Gel Extraction Kit (#K0513).</li> </ul> <p><b>2.4. Incomplete DNA digestion</b> (see 1).</p> <p><b>2.5. Gel shift</b> (see 3.1).</p> <p><b>2.6. Unexpected recognition sites in template DNA.</b> Newly generated target sites in constructed DNA may be overlooked. Re-check your DNA sequence and cloning strategy. Refer to Tables 1.21, 1.22, 1.23 or 1.24 for Newly Generated Recognition Sequences (pp.182-206) to identify all the cleavage sites present in the substrate DNA.</p>
3. Diffused DNA bands	<p><b>3.1. Gel shift.</b> Enzyme that remains bound to the substrate DNA will affect the electrophoretic mobility of the digestion products. FastDigest® Restriction Enzymes BspMI (BveI), HgaI (CseI), AclI (Eco57I), FokI, MboII, TspRI (TscAI) and conventional restriction enzymes AarI, AclI, BclII, BseXI, BveI, CseI, Eco57I, Eco57MI, EcoRII, FagI, GsuI, Lsp1109I, LweI, MboII, MnlI, SclI, TspAI, TscAI, TsoI and TstI are particularly prone to remaining bound to the substrate DNA. This will result in a band or smear above the expected band (see picture below). Heat the digested DNA for 10 min. at 65°C in the presence of 6X DNA Loading Dye &amp; SDS Solution (#R1151) or 0.2% SDS prior to electrophoresis.</p>  <p><b>M</b> – GeneRuler™ DNA Ladder Mix (#SM0331).  <b>1</b> – 0.5 µg λ DNA prepared for loading with 6X DNA Loading Dye (#R0611).  <b>2</b> – 0.5 µg λ DNA prepared for loading with 6X DNA Loading Dye &amp; SDS Solution (#R1151).  <b>3</b> – 0.5 µg λ DNA digested with TsoI (#ER1991), probe prepared for loading with 6X DNA Loading Dye (#R0611).  <b>4</b> – 0.5 µg λ DNA digested with TsoI, probe prepared for loading with 6X DNA Loading Dye &amp; SDS Solution (#R1151).</p> <p><b>3.2. Contaminated reagents.</b> Any restriction digestion reaction components may become contaminated with nucleases due to improper handling or storage. Nuclease contamination causes DNA degradation, which appears as diffused DNA bands on a gel. Perform four control reactions in order to check for the nuclease contamination:  <b>I</b> – without restriction enzyme, <b>II</b> – with a new vial of buffer, <b>III</b> – without restriction enzyme, with a new vial of buffer, <b>IV</b> – with commercially available water e.g. Water, nuclease-free (#R0581).</p> <ul style="list-style-type: none"> <li>Contaminated sample DNA (diffused bands in all controls). Prepare new DNA sample (re-purify DNA).</li> <li>Contaminated enzyme (diffused bands in controls 2 and 4). The enzyme may become contaminated due to improper handling. Use a new vial of enzyme.</li> <li>Contaminated buffer (diffused bands in controls 1 and 4). Bacterial contamination of the reaction buffer will cause DNA degradation. Use a new vial of buffer. Store all buffers at -20°C.</li> <li>Contamination of both enzyme &amp; buffer (diffused bands in controls 1, 2 and 4). Follow the recommendations given above.</li> <li>Contaminated water (diffused bands in controls 1, 2 and 3). Bacterial or DNase contamination in improperly handled water will cause DNA degradation. Use commercially available nuclease-free molecular biology grade water (e.g. #R0581).</li> </ul>