

High Specificity Real-time PCR



Highest quality & performance



Maxima™ SYBR Green/Fluorescein qPCR Master Mix (2X), #K0241

- **Specificity** – Maxima™ Hot Start *Taq* DNA Polymerase and the optimized buffer eliminate non-specific amplification and formation of primer dimers.
- **Sensitivity** – detects low copy number targets.
- **Fluorescein** – supplemented with fluorescein passive reference dye for use with Bio-Rad iCycler®, MyiQ™ or iQ5™.
- **Wide linear range** – accurate quantification across 9 orders of magnitude.
- **Reproducibility and convenience** – ready-to-use 2X master mix.

Product	Cat. #	Size
Maxima™ SYBR Green/Fluorescein qPCR Master Mix (2X)	K0241	200 reactions
	K0242	1000 reactions

Maxima is Fermentas trademark. PureExtreme is registered Fermentas trademark.

Uniform, specific and reproducible results



Fermentas
Part of Thermo Fisher Scientific

www.thermoscientific.com/fermentas

ISO 9001 | ISO 14001 | ISO 13485

MP19